

John Weiss offers 10% introductory discount on the Tecfen range, to celebrate new UK Distributor status

Friday 21 February, 2020

Harlow, UK – 21st February 2020 - [John Weiss & Son](#), the leading manufacturer and distributor of gold-standard ophthalmic micro-surgical instruments and procedures, is proud to be appointed the Exclusive UK Dealer of the Tecfen Medical [corneal trephine blade](#) portfolio and an Authorised UK Distributor for the Tecfen Medical range of [tissue punches](#).

To celebrate this new partnership, John Weiss & Son is pleased to offer a **10% introductory discount** across the **entire Tecfen product range**. This special price promotion is valid immediately and will be available until 31st March 2020. Please call **01279 456256** and quote **TECFENQ120** to receive the 10% discount.

[Tecfen](#) Medical offers a comprehensive selection of corneal trephine blades and tissue punches, including the; complete Corneal Punch System & Components (used with the Long Trephines series for donor or recipient), the Vacuum and Non-Vacuum Trephine Punch series (for donor) and the Radial Vacuum Trephines series (for recipient).

The Corneal Punch System is supplied sterile and includes a Punch Piston and a Piston Holder & Base Plate, where a separately packed sterile and disposable Teflon corneal block fits. The Corneal Punch Piston is designed to hold the Long Trephine Blades when used by hand as a punch to cut the cornea from the donor. The Punch Piston is available in three size; 5-7mm, 7-9mm and 9-11mm.

Long Trephine Blades (for donor or recipient) can be used by hand or with a hand piece, such as the Tecfen Corneal Punch system, and are compatible with most punch systems. They deliver consistent results and range from 2.5 - 18mm in diameter.

Vacuum and Non-vacuum Trephine Punches (for donor) are available as a complete system, with blades ranging from 6mm - 9.5mm in diameter. This provides flexibility of use in instances where the Surgeon is not concerned about the vacuum when cutting the cornea from the donor. Non-Vacuum Trephine Punches are designed with four posts to ensure proper alignment between the punch base and the trephine blade holder, ensuring a perfect cut, every time.

The Radial Vacuum Trephine (for recipient) is a complete vacuum trephine system, with 16 vacuum chambers to provide stability for a clean and precise cut. The recipient trephines range from 6mm-9mm diameter in size and are supplied sterile and individually packaged.

To take advantage of the 10% introductory discount, call John Weiss & Son now on **01279 456256**, quoting offer code **TECFENQ120**. For more information on the Tecfen range, please visit <https://jweiss.co/tecfen>.

Media:



Related Sectors:

Health :: Manufacturing,
Engineering & Energy :: Medical
& Pharmaceutical ::

Related Keywords:

Tecfen :: Trephines :: Surgery ::
Ophthalmology :: Cornea :: Offer
:: Discount ::

Scan Me:



Company Contact:

—

John Weiss & Son

T. 01279 883858

E. rseymour@johnweiss.com

W. <https://www.haag-streit.com/john-weiss/>

[View Online](#)

Additional Assets:

Newsroom: Visit our Newsroom for all the latest stories:

<https://www.johnweiss.pressat.co.uk>